

# A concordance study of the ArcherDx Reveal ctDNA 28 NGS panel and Biocept's Target-Selector™ mutation assay using ctDNA collected in Biocept's CEE-sure blood collection tube

Jason C. Poole<sup>1</sup>, Brian Kudlow<sup>2</sup>, Jill Stefanelli<sup>2</sup>, Skyler Mishkin<sup>2</sup>, Anh Pham<sup>1</sup>, Jeff Chen<sup>1</sup>, Veena M. Singh<sup>1</sup>, Josh Stahl<sup>2</sup>, Lyle J. Arnold<sup>1</sup>.

AACR 2017

<sup>1</sup>Biocept, Inc., San Diego, CA; <sup>2</sup>ArcherDX, Inc., Boulder, CO

Poster # 762

## Background

- The ability to build complex, targeted libraries, free of chemically induced mutations is of major importance for the detection of circulating tumor DNA (ctDNA) from solid tumors found in liquid biopsies. DNA damage occurring as a result of sample handling, polymerase errors or chemical preservatives can create low abundance, falsely positive signals indistinguishable from cancer-derived mutations. A critical component of this, often overlooked, is the proper collection and preservation of the sample paired with the right isolation and detection methodology.
- We undertook to evaluate and compare the performance of samples collected in Biocept CEE-Sure blood collection tubes by generating NGS libraries using the ArcherDX Reveal ctDNA 28 panel. We utilized remnant patient DNA isolated from plasma collected and processed at our CAP/CLIA certified laboratory in San Diego, CA.

## Study Design

Sample	Biocept Result	Mutation Level
1	KRAS Ex2	high
2	KRAS Ex2	low
3	BRAF V600	high
4	BRAF V600	low
5	EGFR L858	high
6	EGFR L858	low
7	EGFR Del19	high
8	EGFR Del19	low
9	EGFR T790M	high
10	EGFR T790M	low
11	Negative	N/A
12	Negative	N/A

Table 1: circulating DNA from patient blood collected in the patented Biocept CEE-sure was extracted from plasma and used in Target-Selector™ assays specific for EGFR (Del19, L858R, T790M), KRAS or BRAF mutations. In certain cases CTC enumeration and biomarker testing was also performed. Sequencing of the amplified Target-Selector™ product was used to confirm the presence of the mutation. The remnant DNA samples were de-identified and sent blinded to ArcherDx and processed through their Reveal ctDNA workflow. **Each sample contained either a high or low mutational burden for one of the 5 clinically actionable markers listed along with two samples which tested negative for all five.**

## CEE-Sure BCT

Figure 1: CEE-Sure™ BCT is designed to both prevent blood coagulation, and to preserve the cells from lysing due to age or transport conditions. CEE-Sure™ allows for room temperature shipment of whole blood for nucleated cell isolation, and for cfDNA and cfRNA detection.

Manufacturing	ISO 9001
Tube Volume	8.0 ml liquid capacity
Tube Type	Plastic or Glass Plastic: 12 months Glass: 24 months
Shelf Life	Glass: 24 months
Storage Temperature	18°C – 25°C
Shipping Temperature	15°C – 30°C
Stability	Blood: Up to 4 days DNA: Up to 8 days


**Advantages**

- Ships at room temperature
- Patented preservative stabilizes cellular and genetic content
- ACD prevents blood coagulation
- Choice of Plastic or Glass
- Used commercially in Biocept's own CLIA lab for 2+ years for CTC capture and ctDNA applications
- Proprietary solution
- Manufactured in accordance with ISO9001


## Testing and Tools

### Target Selector Test

Figure 3: Target Selector is a targeted hotspot mutation test that enriches for rare mutations in a large excess of WT DNA. The test is specific to small regions of interest, highly sensitive, and validated at 7 mutant copies in a background of 14,000 WT at >98% sensitivity. Target-Selector™ probes block WT amplification but allow mutant DNA amplification. A mismatch (mutation) under the switch region creates a large annealing temperature discrimination between the mutant and WT binding temperatures allow mutant specific amplification.



### Target Selector Workflow



### Archer ctDNA Reveal 28

The ArcherDX Reveal ctDNA test is a 28 gene NGS panel that targets key oncogene activating mutations, drug resistance mutations, and yields full coverage TP53. The Reveal ctDNA assay utilizes Anchored Multiplex PCR (AMP) to enrich, tag, and efficiently capture short ctDNA fragments.




Figure 5: ctDNA inputs, naturally fragmented prior to extraction from plasma, are end-repaired, A-tailed, and tagged by ligation to a molecular barcode-containing adapter. Target enrichment with AMP is achieved by two rounds of nested PCR between region-specific primers (GSP1 and GSP2) and a universal primer (P5 Primer). Fully indexed, target enriched libraries are produced in under 7 hours (left). Archer ctDNA Reveal 28 target list (right).

## Biocept Results

### CTC Platform

Sample #	Accession #	Cancer Type	CK+ CTCs	CK- CTCs	ALK	ROS1	FGFR1
1	16-01431	lung	0	31	Not Det.	Not Det.	Not Det.
2	16-01195	colon	-	-	-	-	-
3	16-01350	melanoma	-	-	-	-	-
4	16-01064	lung	0	21	Not Det.	Not Det.	-
5	16-01513	lung	-	-	-	-	-
6	16-01493	lung	-	-	-	-	-
7	16-01494	lung	-	-	-	-	-
8	16-01512	lung	-	-	-	-	-
9	16-01400	lung	-	-	-	-	-
10	16-01414	lung	-	-	-	-	-
11	16-01498	lung	0	11	-	-	-
12	16-01516	unknown	2	24	Not Det.	Not Det.	-

Table 2: Additional information on patients tested in the study from the Biocept Clinical Laboratory using the target selector tumor cell capture platform. Cytokeratin negative CTCs were found in 4 of the patients including the two where no other genetic alterations were found. Highlighted cells indicated signal was detected. Dash lines mean the test was not ordered.

### ctDNA Platform

Sample #	T790M	T790M %	L858R	L858R %	Del19	Del19 %	KRAS	KRAS %	BRAF	BRAF %
1	0/1912	0	0/1912	0	0/1912	0	493/1912	26%	0/2589	0
2	-	-	-	-	-	-	120/14219	<1%	-	-
3	-	-	-	-	-	-	-	-	690/17890	3.86%
4	5/1532	<1%	29/1256	2.2	34/2006	1.6	0/1194	1.2	12/1512	0.79%
5	0/2625	0	33/2592	1.2	0/2325	0	-	-	-	-
6	0/1231	0	1/1229	<1%	0/1231	0	-	-	-	-
7	5/707	<1%	0/712	0	25/686	3.5	-	-	-	-
8	2/1752	<1%	0/1751	0	3/1751	<1.0	-	-	-	-
9	350/13121	2.67%	0/13475	0	1/13474*	0	-	-	-	-
10	1/1342	<1%	0/1343	0	4/1339	0.30%	-	-	-	-
11	0/1336	0	0/1398	0	0/1398	0	0/1481	0	0/1094	0
12	0/5693	0	0/5693	0	0/5693	0	0/7547	0	0/5530	0


Table 3: Complete testing results for the study patients using Biocept's target selector circulating tumor DNA platform. Tabular numbers indicate copy #/mL of plasma tested for either mutant or WT species of DNA. Blue highlight indicates mutation detected. Dashed lines indicate test was not ordered for this patient.

## Archer Results



### Library Quality

Sample	ng/ul	ul Input	ng Input	nM Yield	AVERAGE_UNIQUE_DNA_AND_AMBIG_START_SITES_PER_GSP2	Avg. Unique Reads per GSP2	UNIQUE_FRAGMENT_TOTAL
16-01431	0.107	50	5.35	12.75	121.43	445.69	218015
16-01195	1.04	30	31.2	31.93	165.97	1626.88	495736
16-01350	1.42	20	42.6	150.28	251.3	3582.94	892171
16-01064	0.195	30	5.85	17.10	132.31	474.55	227398
16-01513	0.342	30	10.26	47.93	169.28	1325.01	414318
16-01493	0.154	30	4.62	6.44	107.35	300.56	115502
16-01494	0.176	30	5.28	7.30	89.18	219.43	104696
16-01512	<0.05	35	NA	1.74	42.68	73.23	62279
16-01400	1.21	30	36.3	128.68	246.84	2466.75	681210
16-01414	0.116	30	3.48	6.08	95.96	233.99	120997
16-01498	0.71	30	21.3	14.04	97.82	319.48	146512
16-01516	0.509	50	25.45	55.01	156.49	1333.93	426587

### Detailed Variant Analysis



### Variant Frequency



## Study Results

Patient	Gene	Mutation Level	Biocept	ArcherDx	Biocept Other	Archer Other
1	KRAS	High	X	X	CK- CTCs,	None
2	KRAS	Low	X	X	None	None
3	BRAF	High	X	X	None	None
4	BRAF	Low	X		CK- CTCs, EGFR T790M, L858R, Del19	PIK3CA E545K, MET
5	EGFR L858R	High	X		None	None
6	EGFR L858R	Low	X		None	KRAS G12R, KIT (D572G, T574A)
7	EGFR Del19	High	X	X	EGFR T790M	EGFR T790M
8	EGFR Del19	Low	X		EGFR T790M	KIT (D572G, E583G)
9	EGFR T790M	High	X	X	EGFR Del19*	EGFR (Del19, A750P)
10	EGFR T790M	Low	X		None	None
11	Negative	N/A	X	X	CK- CTCs	KIT K558E
12	Negative	N/A	X	X	CK+ CTCs, CK- CTCs	None

Table 3: Concordance of the tests in 7 of 12 patients. The Biocept test outperformed on sensitivity for the most common markers, but the ArcherDX panel added important additional markers that could guide patient treatment in certain cases. Overall the tests complement one another. For instance one could test the most common markers from Biocept and if negative follow up with a larger panel to look for less common markers that may guide treatment toward therapeutics in clinical trials.

## Conclusions

- In this pilot evaluation, the tumor status of 12 patients were tested using both the Biocept Target Selector Assay and the Archer ctDNA reveal 28 NGS panel.
- Archer was able to produce high quality libraries from the stored sample, indicating that little DNA damage occurred during preservation in the CEE-Sure tube.
- The two methods were shown to have good concordance and while the target selector tests excelled at finding low abundance markers, the Archer panel uncovered mutations untested in the Biocept profiling.
- This study demonstrates that both breadth of coverage and depth of analysis are important factors for today's oncologist when considering which diagnostics to perform on their patient. The rise of the liquid biopsy provides an powerful, informative, non-invasive option when looking for actionable biomarkers.

